


Blood Pressure Cuff DO's and DON'Ts

DO



✓ CLEAN CUFF IN GOOD CONDITION

DON'T



✗ CUFF WITH FUR COVERING THE HOOK AND LOOP SELF-FASTENER



✓ CUFF PLACED WITH TRANSDUCER TUBING DIRECTED AWAY FROM THE CAT



✗ CUFF PLACED WITH TRANSDUCER TUBING DIRECTED TOWARD THE CAT



✓

- CUFF POSITIONED AT HEART LEVEL
- CUFF PLACED ON NONDEPENDENT LEG IN LATERAL RECUMBENCY
- TRANSDUCER TUBING AWAY FROM THE CAT



✗

- CUFF POSITIONED BELOW HEART LEVEL
- CUFF PLACED ON DEPENDENT LEG IN LATERAL RECUMBENCY
- TRANSDUCER TUBING TOWARD THE CAT

For more information, visit catvets.com/hypertension-toolkit.

© 2021 American Association of Feline Practitioners

